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
Roughness measurement is of critical importance in implantology since it directly affects the acceptance by the human body. For example here a measure of a dental implant surface roughness at the root of a thread.

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Size, shape, and surface finish...

The quality of measures in the departments of metrology and laboratories is constantly drawn upwards, furthermore new standard of roughness measurement taking surface into account and no longer just profile involves changes in the means of measurement. We discovered InfiniteFocus, a device capable of measuring automatically the three elements of dimensional precision, shape accuracy and quality of surface finish. Meeting with Mr. Jérôme Orteu, Director of Microcube and exclusive reseller for the French-speaking Switzerland of this extraordinary device produced by the Austrian company Alicona.

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

In the field of Microtechnology, miniaturization and combination of functions on parts are always more topical, what compels companies to create more and more complex parts, also more accurate with tighter tolerances. Measurement tools have to follow this trend and for companies, changes at this level can be of heavy consequences, not just in budgetary terms. With its InfiniteFocus system Alicona responds to this need in a very interesting way. A single investment in this tool replaces the purchase of two complete measuring systems (one for shapes and sizes and the other for surface finishes).

Ergonomic equipment

The measurement device includes a gantry which supports the lenses and lighting means (Z axis), a motorized cross-table (X/Y axes), gripping devices of parts, screen, control and a cabinet including the computer, the engines and control systems. According to the needs, additional options can be set up on the table, for example a device for automatic rotation of parts to gather 3D data on 360° (to compare with a 3D body). The whole design was thought "operator" and working on the system is very easy (ask for a demonstration, the operation is really very intuitive). The scan is performed once, then the operator can work with all data.

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A simple and proven technique

The principle of operation of the machine is based on an optical measurement using the principle of focus variation. The part is positioned on its support, the operator controls the three axes with a joystick. In the case he positions approximately the optics in front of a revolution part raised vertically, the device recognizes the axes and automatically focuses on the central hole (for example). Then the operator must indicate the extreme measuring points (the lowest and highest points), this simply using visual focuses on the screen of the computer by varying the height of the focus. Once two points registered, scanning is automatic and the screen displays a 3D Image of the part. Tabs contain all the data. For the scanning, the lenses moves vertically and the device records all information about the different layers of clear focuses.

Various lenses, same strengths

The system is delivered with a standard rotating head comprising three lenses to choose between 2.5 x and 100 x. In all cases, the technology offers the considerable advantage of being able to measure while keeping the front lens well away from the part. With a 5 x lens, the measuring distance is 23 mm, which means that one can measure relatively deep interiors. The measurement area with a 2.5 x lens reaches 5,7 mm x 4,3 mm. With a 100 x lens, it reaches 150 µm x 100 µm. The resolution has to be chosen based on the aims of the measure, of course the higher it is, the more scanning takes time. Resolution is usually set between 200 and 400 nm; for a roughness measure it can be lowered to 10 nm!

Powerful automation

Often the parts to measure are larger than the above mentioned areas. In this case, the operation is simple, one indicates the extreme points to measure (x and Y and always Z-, and Z +) and the system automatically performs scanning path to cover all of the area to be measured. If sections are of no importance, it is also possible to tell the system not to measure on these areas. In the case of measures of batches of parts, the system will automatically pass from one to the next (for example when using a base comprising many parts on the cross-table) and once measurement is complete, the operator has instantly access to individual and global reports. Many embedded algorithms provide all the statistical information imaginable.

The different functions of the InfiniteFocus device

- Measurement of profiles roughness
- Measurement of surfaces textures
- Measurement of 2D shapes
- Measurement of volumes
- 2D Measures
- Measurement of 3D shapes
- Measurement of contours
- Measurement of differences
- Measurement of cutting tools

True color information

In the case of a broad surface being scanned, the repositioning of the different scanned areas is, of course, based on the geometrical features of the part and the moves of axes, but also using the color information of pixels, thus alignment is perfect and there is no risk of deviation due to a "bad overlay". At any time the operator can modify the light, he has the possibility to combine a coaxial source of white light and a ring of LEDs around the lenses. It can also add a polarizing filter for best results depending on the material to measure.

A few facts microcube

Foundation: 2006

The purpose of the company: service provider in the field of industrial automation, measurement and software development.

References: Federal Office of Metrology, large watchmakers, pharmaceutical and food groups, mainly in Switzerland.

Strong points: ergonomics and quality of user interfaces, versatility and speed, marriage of metrology and automation skills

Total security for clients (all source codes are given to them at the end of a project)

Alicona certified service partner (www.alicon.com), exclusive distributor for the French-speaking Switzerland.

All types of parts

The maximum measurement surface is 100 mm x 100 mm (1000 x 1000 with InfiniteFocus-Xlarge) and vertical move of the Z axis is also 100 mm. According to the needs, the height of the gantry is adjustable in order to position a very high piece on which the measure should be carried out. The principle of operation remains always the same. This change in height is normally made ex-factory but in the case of a fundamental change of parts to measure in a company, a retrofit operation is possible by Mr. Orteu in French-speaking Switzerland and by Alicona's specialists in other countries. The measurement system is also able of processing parts showing very small radiuses and high slopes (up to 85 °).



Measurement traceability

To ensure the complete traceability of measurement, Alicona has worked with the PTB (Physikalisch-Technische Bundesanstalt) that is the National Agency for Metrology in

Germany. Calibration tools are certified by the PTB and used for the calibration of the system. It works with a statistically known level of uncertainty which guarantees the accuracy of measurements.

And this new standard of roughness?

ISO 25178 standard defines the quality of 3D surface finishes. Within a short time span, the well known Ra states are likely to be replaced by the new Sa states. This new standard (and new way of seeing things) allow to take into account the new means of measuring without contact. Alicona specialists participated in the development of this standard, InfiniteFocus is a tool that really corresponds to the needs of this specific field.

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The whole design was thought "operator" and working on the system is very easy.

InfiniteFocus, whom for?

Aside from the benefits developed in this article (simplicity of operation, "two in one" device, diversity and power of measures or even perfect integration in the evolutions of the ISO standards), this device offers a combination of unique performances of interest to all areas where the "small and accurate" is requested. Namely microtechnology, microfluidics, microelectronics, plastic injection, tribology... and even forensic medicine. Mr. Orteu says to conclude: *"With InfiniteFocus, we are finally able to quantify certain elements so far completely dependent on human eyes. As an example we are setting up systems to validate visual quality in top watchmaking that will be able to replace the current visual controls. The system is capable of guaranteeing quality, even when it's merely aesthetics".*

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cl055841GE © Laurent Bouvet / Agence Rapsodia
Alicona Imaging GmbH